

HOW TO USE A

Dry Powder Inhaler — Diskus[®]



A Diskus[®] is an inhaler shaped like a hockey puck or flying saucer. It delivers pre-set doses of medicine in powder form. The medicine gets to your airways when you take a deep, fast breath in from the inhaler. To keep your asthma under control, it is important to take your medicine as prescribed by your doctor or other health care professional and to use the proper technique to deliver the medicine to your lungs. If you don't use your inhaler correctly, you won't get the medicine you need.

Here are general steps for how to use and clean a Diskus[®] inhaler. Be sure to read the instructions that come with *your* inhaler. Ask your doctor, pharmacist, or other health care professional (such as nurse practitioner, physician assistant, nurse, respiratory therapist, or asthma educator) to show you how to use your Diskus[®] inhaler. Review your technique at each follow-up visit.

1. Holding the Diskus[®] flat (like a flying saucer) in one hand, place the thumb of your other hand in the thumbgrip and slide the cover open until it clicks.
2. With the mouthpiece facing you, slide the lever away from your until it clicks. This will load a dose of medicine. Do **not** tip or shake the inhaler.
3. Stand up or sit up straight.
4. Take a deep breath in and blow out completely to empty your lungs. Do **not** blow into the inhaler.
5. Place the mouthpiece of the inhaler in your mouth and close your lips around it to form a tight seal.
6. Take a fast, deep, forceful breath in through your mouth.
7. Hold your breath and count to 10.
8. Take the inhaler out of your mouth. Breathe out slowly, facing away from the inhaler.
9. Slide the cover closed.
10. If you are supposed to take more than 1 inhalation of medicine per dose, wait 1 minute and repeat steps 1 through 9.
11. If using an inhaled corticosteroid, rinse out your mouth with water and spit it out. Rinsing helps to prevent an infection in the mouth.
12. Store the Diskus[®] in a cool, dry place (not in the bathroom).

Note: A Diskus[®] inhaler has a built-in counter to tell you how many doses are left. When the counter gets to "0," throw it away. Be sure to get a refill of your inhaler before the counter gets to "0."

HOW TO CLEAN A DISKUS[®] INHALER

- Wipe the mouthpiece at least once a week with a dry cloth.
- **Do NOT use water** to clean a Diskus[®] inhaler.



National Asthma Education
and Prevention Program (NAEPP)

Convened by the
National Institutes of Health